


Remote control presentation



www.chromatotec.com

Remote control advantages

- Display of analyzer window
 - Allows full control of Chromatotec™ softwares
- Possibility to give us the hand on the analyzer if you have any problem
- Very useful for ATEX zone

Team Viewer

- What is it?
 - Software for remote control
 - Similar to UltraVNC
 - Allows control of analyzer all over the world

Team Viewer VS UltraVNC

- Less configuration
 - Possible to use it with the majority of firewall and proxy servers
- Better safety
 - In compliance with actual norms
 - Not possible for an intermediary server to see the data during connection




Team Viewer VS UltraVNC

- Possibility to connect even with very slow environments
- Possibility to record video of your work (training)
- Easier data transfer

Team Viewer VS UltraVNC

- Windows resolution adjustable



Team Viewer VS UltraVNC

- Temporary or fixed password



- Multilanguage software


Team Viewer VS UltraVNC

- Compatibility:
 - Windows 8 / 7 / Vista / XP / 2000
 - Windows Server 2012 / 2008R2 / 2008 / 2003 / 2000, Windows Home Server/Home Server 2011
 - Mac OS X 10.5 and later (including Mountain Lion)
 - Linux: Red Hat, Fedora, Suse, Mandriva, Debian and Ubuntu
 - iPhone/iPad
 - Android

Team Viewer Licence

- Team Viewer Host (on analyzer) is free
- Team Viewer (which allows you to control another computer) **needs a licence**
 - With 1 licence, you can control up to 3 Team Viewer Hosts
 - If you want more, you need to pay an additional cost

Team Viewer principle



Thanks a lot for your attention



www.chromatotec.com

11